

BABY DRAPE




SURGIWEAR has developed this drape after feedback from surgeons. It offers many advantages:

- Drape is sterilized in a CE certified facility.
- No transfer of fabric fibers into the wound.

Use of this drape will result to:

- No or less post operative swelling of wound.
- Negligible post operative fever.



- Supple wound area.
- Minimum post operative pain.
- Lesser incidence of post operative infections.
- Lesser morbidity and faster recovery of patient.
- **Ultimate cost saving by way of lesser use of medicines.**

Uses: For all kind of surgery on baby, it helps in maintaining body temperature. It is used in many places as low cost drape

d D201

Material: Translucent soft poly film.

Sheet Size 1.2 x 1.2 m

Adhesive area 25 x 20 cm with soft adhesive

Supplied ready to use in peel open pouches.
Sterilized by EO

be all free from disease & sufferings

SURGIWEAR®
AN ISO 13485 : 2003 COMPANY