

# GOWNS

SURGIWEAR has a range of comfortable gowns. The range includes full gowns, half gowns and gowns with level 4 protection.

- Gowns are sterilized in a CE certified facility.
- Minimum transfer of fabric fibers into the wound.
- No chemicals present in the fabric of the gown.

Use of these gowns will result to:

- No or less post operative swelling of wound.
- Negligible post operative fever.
- Supple wound area.
- Minimum post operative pain.
- Lesser incidence of post operative infections.
- Lesser morbidity and faster recovery of patient.
- **Ultimate cost saving by way of lesser use of medicines.**



Uses: For wearing during surgery

## *h* HALF GOWN D800

**Material:** Plain polyethylene sheet 25 micron thick

## *h* HALF GOWN D810

**Material :** White non woven SMS 25 gsm water repellent.



## *f* FULL GOWN D802

Wrap-around Gown. Level 2 protection as per ANSI / AAMI Standard

**Material:** Water Repellent Spun Lace non woven 65-77 gsm

## *f* FULL GOWN ECO E812

Wrap-around Gown. Level 1 protection as per ANSI / AAMI Standard

**Material:** Water Repellent Blue non woven SMS 40-50 gsm

## *f* FULL GOWN BVB D805

Wrap around gown critical areas made of Ahlstorm's Breathable Viral Barrier fabric. Provides level 4 protection in critical areas as per ANSI and AAMI standards. Protection level certified by Ahlstorm.

Front and part of sleeves made of BVB Ahlstorm fabric and rest made of 66-70 gsm water repellent spun lace fabric

be all free from disease & sufferings

**SURGIWEAR®**  
AN ISO 13485 : 2003 COMPANY