

# NEUROSURGICAL PATTIES

## *i* INTRODUCTION

Patties with thread are used to mop blood, fluids and as protective cover for delicate tissues in Neurosurgery and many other specialties. Purpose of thread is to identify the inserted patty.


SURGIWEAR has developed world class patties. These are made of spun laced / thermally bonded non woven fabrics. Most of them have viscose as main fiber. These are highly absorbent.

## *a* ADVANTAGES

- Mechanically/thermally bonded fibers reduce shedding of fibers, reduced linting.
- Bonded fibers make removal of patties easy.
- Radio-opaque types can be identified in x-ray pictures.

## *d* DESCRIPTION

Non radio-opaque types have black nylon thread bonded with the patties. Radio-opaque types have 'blue polypropylene with bariun sulphate' thread bonded with the patties. SURGIWEAR has developed a special process to bond the thread with patties. The Special automatic machine, developed in house by SURGIWEAR, can bond and cut many lines of patties automatically.



Currently three types of base material and two types of thread are used. Thus there are six types of patties.

Material	Prefix	Example
Thin	No prefix	WM—
Thick	T Prefix	TWM—
Soft thin	P Prefix	PWM—
Thread	Suffix	Example
Black Nylon	No	-WM—
Blue PP R-opaque	R	-WM—R

## *p* PRESENTATION

Neurosurgical patties are supplied in sterile peel open packs on a holding card for easy handling. Following sizes are regularly available, other sizes can be manufactured on demand.

Code no.	Size (WxL)	Number in a pack
WM0505	5 x 5 mm	10 patties
WM1010	10x10 mm	10 patties
WM1015	10x15 mm	10 patties
WM2020	20x20 mm	10 patties
WM2030	20x30 mm	10 patties

be all free from disease & sufferings

**SURGIWEAR®**  
AN ISO 13485 : 2003 COMPANY

WM2040	20x40 mm	10 patties
WM2050	20x50 mm	10 patties
WM5070	50x70 mm	3 patties
WM70100	70x100 mm	2 patties

### **WARNING**

The Thread is for identification only.

### **PRODUCT INFORMATION DISCLAIMER**

G. Surgiwear Limited has exercised reasonable care in the choice of materials and manufacture of this product. G. Surgiwear Limited excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for a particular purpose. G. Surgiwear Limited shall not be liable for any incidental or consequential loss, damage or expense, directly or indirectly arising from use of this product. G. Surgiwear Limited neither assumes nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product.

**SURGIWEAR®**  
AN ISO 13485 : 2003 COMPANY

European Authorised Representative  
Obelis s.a.,  
34, Av. De Tervuren, bte 44 BELGIUM  
Tel: +32 2 732 5954 Fax: +32 2 732 6003  
e-mail: mail@obelis.net

Made in India by:  
**G. Surgiwear Limited**  
Rasoolpur Jahanganj, Shahjahanpur-242001  
Tel: +91 5842-223818 Fax: +91 5842-222190  
e.mail: surgewear@hotmail.com  
Delhi: 65454565, Mumbai: 66994258, Ahmedabad: 27545218

[www.surgiwear.co.in](http://www.surgiwear.co.in)