

# O-SCOPE DRAPE

## *i* INTRODUCTION

Use of operating microscope is increasing every day. Operating microscopes can not be sterilized. Formaldehyde tablet system can reduce some bacteria, but can not sterilize an equipment. The best answer to this problem is use of a good drape made specifically for the op-microscope.

## *a* ADVANTAGE

Surgiwear has designed microscope drape in the name of O-Scope drape for specialist surgeons. All sides of microscope including lens is covered.

## *S* SPECIAL FEATURE

Special features of this drape are following:

The lens cover is at an angle. This avoids direct reflection of light into the observing eye.

The lens cover is made of special plastic sheet with high clarity.

It is suitable for all kinds of microscopes.



## *d* D204

Size 1.25 x 2.5 mt

Plain polyethylene sheet tube 30-40 gsm, one objective lens cover & three pair eye lens cover with two tying straps.

**Uses:** For use on operating microscope.

be all free from disease & sufferings

**SURGIWEAR®**  
AN ISO 13485 : 2003 COMPANY