

# WHITE MOP



Cotton Mops are used in all kinds of surgical operations. Conventionally mops are made of surgical cotton gauge. The biggest disadvantages with these mops, is that these carry chemicals used during manufacturing process. Moreover recycled mops carry protein debris of last operation and all the chemicals used by washer man during washing.

SURGIWEAR has developed WHITE MOP surgical sponges. Specially designed for surgeons, for almost every kind of surgery. It will keep your patient's wound free from chemicals.

## *a* ADVANTAGES

- WHITE MOP is chemical free. No unwanted chemicals will go into the wound.
- WHITE MOP is less linting during surgery. Their use will minimize fiber granulomas, post operative infections due to foreign body reaction to shredded fibers.
- WHITE MOP is radio opaque. During the manufacturing process, a radio opaque thread containing Barium sulphate is added to make it radio-opaque. This quality will help surgeon to identify any lost sponges.
- WHITE MOP is super absorbent. The rate of absorption is very fast.
- WHITE MOP is meant for single use only.

## *p* PRECAUTIONS DURING USE

WHITE MOP should be moistened before use. Avoid rubbing it on rough surface.

**Presentation:** Supplied in peel open pouches, ready to use & sterilized by EO

**Material :** Non woven spun lace perforated 100 %, cotton of approx 40 gsm

Product Name	Code	Size
White Mop Extra Large	WM101	5 x 30 cm
White Mop Large	WM102	30 x 20 cm
White Mop Extra Large	WM103	25 x 17 cm
White Mop Medium	WM104	20 x 15 cm
White Mop Small	WM105	15 x 12 cm
White Mop Patch	WM106	12 x 10 cm
White Mop Stripe	WM107	4 x 15 cm
General Surgery Set Large	WM001	
General Surgery Set Medium	WM002	
General Surgery Set Small	WM003	
Neuro Surgery Set Large	WM011	
Neuro Surgery Set Small	WM012	

be all free from disease & sufferings

**SURGIWEAR®**  
AN ISO 13485 : 2003 COMPANY