

PLAIN DRAPE

SURGIWEAR has developed these drape after feedback from surgeons. It offers many advantages:

- No fluid can pass through the drape in operating area.
- Drape is sterilized in a CE certified facility.
- Minimum transfer of fabric fibers into the wound.
- No chemicals present in the fabric of the drape.

Use of this drape will result to:


- No or less post operative swelling of wound.
- Negligible post operative fever.
- Supple wound area.
- Minimum post operative pain.
- Lesser incidence of post operative infections.
- Lesser morbidity and faster recovery of patient.
- Ultimate cost saving by way of lesser use of medicines.

Uses: For Covering patient body & instrument trolleys during surgery

| Product Name | Product Code | Size |
|--------------|--------------|------|
|--------------|--------------|------|

Material: 25 micron thick PE hygiene film

| | | |
|-------------------|------|------------|
| Plain Sheet Small | D300 | 150X120 cm |
| Plain Sheet | D301 | 210X120 cm |



| | | |
|-------------|------|-----------|
| Plain Towel | D303 | 66X120 cm |
|-------------|------|-----------|

| | | |
|----------------------------|------|-----------|
| Mayo's Trolley Cover Small | D304 | 120X66 cm |
|----------------------------|------|-----------|

| | | |
|----------------------------|------|------------|
| Mayo's Trolley Cover Large | D305 | 110X100 cm |
|----------------------------|------|------------|

Material: SMS Polypropylene Water repellent 50 gsm

| | | |
|----------------------|------|------------|
| Plain Sheet Small NW | E310 | 150X120 cm |
|----------------------|------|------------|

| | | |
|----------------|------|------------|
| Plain Sheet NW | E311 | 210X150 cm |
|----------------|------|------------|

| | | |
|----------------|------|-----------|
| Plain Towel NW | E313 | 120X75 cm |
|----------------|------|-----------|

| | | |
|-------------------------------|------|-----------|
| Mayo's Trolley Cover Small NW | E314 | 130X75 cm |
|-------------------------------|------|-----------|

| | | |
|-------------------------------|------|------------|
| Mayo's Trolley Cover Large NW | E315 | 150X123 cm |
|-------------------------------|------|------------|

Material: Spun lace Non woven

| | | |
|-----------------|------|------------|
| Neo natal sheet | D504 | 75 X 80 cm |
|-----------------|------|------------|

Material: Spun lace Non woven laminated with PE 90 gsm

| | | |
|-----------------|------|------------|
| Plain Sheet Int | I321 | 210X120 cm |
|-----------------|------|------------|

| | | |
|-----------------------|------|-----------|
| Plain Sheet Small Int | I323 | 120X75 cm |
|-----------------------|------|-----------|

| | | |
|-----------------|------|----------|
| Plain Towel Int | I327 | 60X75 cm |
|-----------------|------|----------|

be all free from disease & sufferings

SURGIWEAR®
AN ISO 13485 : 2003 COMPANY